

# Auto Logic 110 Mattress Overlay

ARJOHUNTLEIGH  
GETINGE GROUP



Auto Logic® 110 mattress overlay system is part of the logic range from ArjoHuntleigh. *Auto Logic 110* is recommended for round-the-clock prevention and management of up to grade 4 pressure ulcers with frequent monitoring, nursing intervention and repositioning.

Forget about adjusting pump settings during installation. For optimal pressure area management of very high risk patients, the *Auto Logic* pump uses Self Set Technology (SST) to control the mattress pressure.

In both Active (Alternating) and Reactive (Constant Lower Pressure) therapy modes, an advanced microprocessor system regularly assesses the body mass distribution of patients and readjusts cell air pressures to suit their individual requirements.



The **dual-therapy pump** provides leading pressure management therapy options:

- Active (Alternating) mode
- Reactive (Constant Lower Pressure) mode

The **AutoFirm** function provides a temporary stable support surface to assist patient transfer and other nursing procedures.

...with people in mind

# Key features



**1** CPR Control

- Rapid deflation to assist resuscitation
- Simple one handed operation
- Positioned at head end



**2** Cell in Cell (CIC) System

- Centralises patient on the mattress
- Increases stability of edge to aid transfer off and on the surface



**3** Battery option

- Allows continuity of chosen therapy for a minimum of 8 hours during patient transport or in the event of a power cut



**4** Transport Facility

- Simple operation via disconnection or tube-set from the pump
- Provides patient support for up to 12 hours



**5** Micro Air Loss

- Helps reduce humidity at the patient - mattress interface



**6** Cover Options

- Standard Dartex® Cover
- Advantex® Cover - lower friction properties
- eVENT® Fabric Cover - lower friction and moisture vapour permeability



**ARJOHUNTLEIGH**  
GETINGE GROUP

Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies.  
© ArjoHuntleigh, 2015

■ ArjoHuntleigh AB, Hans Michelsensgatan 10, 211 20 Malmö, Sweden  
Phone: +46 (0) 10 335 45 00

[www.ArjoHuntleigh.com](http://www.ArjoHuntleigh.com)

**GETINGE GROUP**

**Getinge Group** is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of **ArjoHuntleigh**, **Getinge** and **Maquet**. **ArjoHuntleigh** focuses on patient mobility and wound management solutions. **Getinge** provides solutions for infection control within healthcare and contamination prevention within life sciences. **Maquet** specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.